



## Certificate of analysis

<b>Product Description</b>	Goat anti-Human IgE HRP Conjugate
<b>Cat. No.</b>	M10198798
<b>Lot No.</b>	D8397
<b>Host species</b>	Goat
<b>Concentration</b>	0.25 mg/ml
<b>pH</b>	7.2
<b>Form</b>	In 0.05 M KPO4 buffer, 0.85% NaCl, containing proteins, preservatives and 50% Glycerol.
<b>Storage</b>	-20°C, avoiding multiple freeze thaw cycles
<b>Dilution Factor</b>	Greater than 1:1,000 in the human IgE ELISA assay.

Important Note: This product as supplied is intended for research use only. not for use in human, therapeutic or diagnostic applications.