



CERTIFICATE OF ANALYSIS

Normal Monkey Serum

Product Description	Pooled normal rhesus Monkey Serum
Product No.	M1035206
Description	Normal rhesus monkey serum is used to study the normal protein composition. monkey serum is used for blocking in antibody-based applications or for use in cell culture media.
Lot No.	J0524C
Format	Liquid
Species of Origin	rhesus Monkey
Preservative	None
Total protein	5.8g/dL
Hemoglobin	normal
IgG	normal
pH	7.4
Sterile	no
Storage	Store at 2-8°C for one week and -20°C for long time. Avoid multiple freeze/thaw cycles.
Expiration Date	June 2029

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.