

Certificate of analysis

CRP Free Serum

| | |
|--|--|
| Product No. | MK1015155 |
| Description | Liquid, Stabilized, CRP Deficient Serum. |
| Lot | 35649 |
| CRP Concentration | 0.03 mg/L |
| Format | Liquid |
| Preservative | 0.1% Sodium Azide |
| Total protein | 6.4g/dL |
| Filtration | 0.2 μ m |
| Application | CRP free serum can be used for preparation of CRP standards and calibrators. |
| Date of expiration and/or retest date | February 2029 |
| Storage | \leq -20°C, Avoid repeated freezing and thawing. |

For research and further manufacturer use only,
Not for Human or Animal therapeutic