

Certificate of analysis

CRP Free Serum

Product No.	M42155
Description	Liquid, Stabilized, CRP Deficient Serum.
Lot	35649
CRP Concentration	0.03 mg/L
Format	Liquid
Preservative	0.1% Sodium Azide
Total protein	6.4g/dL
Filtration	0.2 μ m
Application	CRP free serum can be used for preparation of CRP standards and calibrators.
Date of expiration and/or retest date	February 2029
Storage	\leq -20°C, Avoid repeated freezing and thawing.

For research and further manufacturer use only,
Not for Human or Animal therapeutic